1. Allgemeine Informationen

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1. General Information

Product description and package contents

The CCG® system consists of CCG® and CCG®-GF bands as well as CCG® stabilisers. CCG bands and CCG® stabilisers are patented products. The CCG® band is formed into a loop and compressed using a special cerclator. The CCG® stabiliser is a curved plate with dentiform protrusions on both sides which penetrate into the bone, thus providing a movement-free bone support. CCG® stabilisers must only be used together with CCG® bands tied around them. CCG® bands and CCG® stabilisers are made of pure titanium according to ISO 5832-2 and ASTM F 67. The titanium surface is roughened, enabling osseointegration. The batch or lot number(s) of the implants used must be documented in the patient records. Labels are included in the packaging of the implants for this purpose.

Information on reuse

A previously used band which has had to be removed must not be reused because of the deformation which has occurred. In this case, the cerclage must be repeated using a new band.

Warnings and precautions

These products may only be used in their original condition or with modifications expressly described in the surgical technique. For the successful implantation of CCG® compression bands and CCG® stabilisers it is of course a requirement that the surgeon be fully conversant with the methods of osteosynthesis when treating bone fractures and the problems associated with revision operations of both cemented and cementless hip arthroplasties.

In general the implants are sterilized by means of gamma rays (at least 25 kGy). Implants supplied "sterile" may not be re-sterilized by the buyer.

Key to the symbols used on the product label

	Do not reuse
	Do not re-sterilize
	Lot number
	Sterilization by irradiation
	Quantity
	Caution
	Order number
	Expiry date
	Manufacturer

2. Storage and handling

Storage

Implants must always be stored unopened in their protective packaging.

The implants are packed in such a way that they can be stored at normal room temperature and relative humidity (corresponding to normal storage and working environment in our latitudes; this means about +15°C and +35°C, depending on the time of year, and corresponding normal relative humidity) without risking any impairment of the packaging, the sterility, the product, etc.

Do not expose sterile implants to sunlight without protection (i.e. without protection from light by the packaging carton).

Handling

CCG® Bands and Stabiliser

Before using the sterile implant check the package for damage, which may have compromised product sterility.

If the original packaging is damaged the implant is no longer sterile and must not be implanted. When unpacking the implant, check that it corresponds with the description on the product label (part number and size).

When removing the implant from the packaging, appropriate aseptic procedures must be observed by the surgical staff. The outer and, if applicable, the middle PE pouch/blister can be removed by the non-aseptic surgical assistant; the innermost pouch/blister must be removed by

the aseptic surgical assistant. The implant must not come into contact with any objects that might damage its surface. Before insertion, every implant must be visually checked for damage.

Use of the implants of the CCG® system is permitted only with the corresponding instruments.

3. Instructions for use - special notes an use

Patient information

Prior to surgery, every patient must be informed in detail about the operation and possible adverse effects.

A. Surgical technique CCG® bands

1. Reposition the bone fragments to be fixed using reposition forceps. Repositioning must not be carried out with the CCG® band.
2. Using a curved bone rasp, expose the bone superiostally at the cerclage position.
3. Insert the CCG® loopawl from above so that it maintains close contact with the bone. If using the CCG®-GF: Slide the moveable spike into the appropriate position.
4. Hook in the CCG® band so that the clasp is on the outside and pull through.
5. Pull CCG® compression band through the clasp and tighten.
6. Bring the carriage of the CCG® cerclage tightener all the way to the front by turning to the left, then tighten the band and fixate it with the lever.
7. The compression is carefully carried out by turning to the right with thumb and index finger only. The intended compression has been achieved when the band cannot be moved.

Caution: Using excessive force can result in undesirable over-compression! When the gap of the clasp begins to open, the permissible compression has already been slightly exceeded. The compression must not be continued at this point.

8. Bend the cerclage tightener with the band slowly up to 90°.

WARNING

During bending with the cerclage tightener, a combination of tensile and bending forces affect the CCG® band. Therefore the bending may only be performed up to 90° so that the band will not slip back. Quick bending over 90° will increase the tensile and bending forces, which increases the strain, and the CCG® band may break. Once the CCG® cerclage tightener has been removed, the bending is completed by the cerclage fixator according to the instructions, using only bending, but not tensile forces. This eliminates the risk of breakage. This handling instruction is very important and must be followed.

9. Loosen the lever and remove the cerclage tightener.
10. Insert the CCG® fixator into the CCG® band, bend the band, and shorten it. Use of the CCG® fixators eliminates the risk of band breakage during bending, and it prevents an incorrect shortening of the band. The wire cutter is not a substitute for the CCG® fixator.
11. Set the fixator with one end onto the end of the band and impact it firmly. For application with wire cutter and pestle: Instead of the CCG® fixator, the CCG® band is cut off by a wire cutter 1 cm above the clasp. Then the protruding end of the band is bent by a pestle and impacted onto the clasp.

B. Surgical technique CCG® stabiliser

1. Only use the CCG® stabiliser in combination with CCG® cerclage bands.
2. Choose the correct length of stabiliser according to the cortical bone to be supported or the shape of the fracture line.
3. The stabiliser is applied longitudinally, axially to the bone, and secured using reduction forceps. Bone reduction should be carried out using the reduction forceps, and not by the compressive effect of the CCG® bands.
4. The number of stabilisers used will depend on the cortical defect to be bridged. The positioning and spacing of the stabilisers will have to be adapted to the intra-operative situation.
5. When using several CCG® stabilisers, these should be placed parallel to one another. Contact between the stabilisers must be avoided in all cases.
6. Furthermore, the stabiliser should be positioned in such a way that direct contact with other metal implants is excluded.
7. Use a tappet to tap the stabilisers protrusions evenly into the bone.
8. Secure the stabiliser into place with CCG® bands.
9. When using two CCG® cerclage bands, position one band proximally, the other distally to the reposition forceps. When using more than two CCG® cerclage bands or longer stabilisers, the use of an additional reposition forceps is recommended.
10. The CCG® cerclage bands must be placed perpendicularly to the stabiliser between two grooves to ensure optimum connection between CCG® cerclage band, stabiliser and bone.

Indication

The CCG® system is suitable for particular forms of osteosynthesis, particularly for cerclage of the greater trochanter and femur associated with revision to hip arthroplasties. In contrast to cerclage wires, the CCG® compression band allows a functional compression to be achieved. It is

assumed that an implant revision is already indicated in this connection.

The CCG® system fulfils a supporting function. The following indications apply:

- Reattachment of the greater trochanter
- Bone support in case of a damaged proximal femur (e.g. by bone cement, loosened stem, or polyethylene wear particles), with a possible cancellous bone graft in the proximal femur area.
- Cerclage of bone fissures.
- Fracture protection during rasping of a new prosthetic canal and prior to implantation of a new stem.
- Closing of cortical fenestration.
- Fractures in the area of the femoral prosthesis stem.
- To improve bone-implant contact.

In general, CCG® and CCG®-GF have the same function. In the case of conical shape of the bone, primarily above the lesser trochanter, the CCG®-GF is indicated rather than the CCG®. CCG® stabilisers are specially indicated in fractures and shaft fissures in which the use of CCG® bands alone would not grant sufficient stability.

Contraindications

An implant of suitable design and length is a prerequisite for a successful use of the CCG® system in hip revisions. Complications such as unstable implantations, implant subsidence, stem and prosthesis breakage cannot be prevented entirely by the CCG® system. Therefore CCG® is not indicated in such cases. The CCG® system is not indicated in cases with large defects of the cortical bone and hence unavoidable contact with the implant which would lead to metal wear.

Complications

CCG® band may break as a result of improper handling. Please read the instructions. Complications which result from incorrect indications, surgical technique or aseptic conditions remain the responsibility of the operating surgeon and cannot be attributed to either the manufacturer or the supplier of CCG®.

Information about the Risks of Interfering Interaction

Undesired effects that harm the patient can occur when performing MRI / CT scans. Possible effects include artefacts, warming of the implant, induction of electrical currents, loosening of the implant. The manufacturer's instructions for use must be studied before the device is used. In case of doubt, reference implants must be tested for suitability in the respective MRI / CT instrument within the scope of an individual risk assessment. The patient must be informed about the risks.

Avoidance/minimization of side effects

A tolerance test should be performed before the operation in the case of any uncertainties about the patient's tolerance of the implant material. The appropriate allergy test patches and test instructions are available from your local sales office.

Postoperative treatment

In accordance with the surgical technique. Use recognized procedures and techniques in postoperative care and treatment.

- Document postoperative treatment in accordance with the hospital's standard procedures, e.g.:
- surgical records
- postoperative X-rays
- patient's card
- regular checks and follow-up examinations.

4. Cleaning, Disinfection and Sterilisation / Resterilisation of Instruments

For information about the treatment of the instruments see the Cleaning and Treatment Instructions Lit. No. 860501 or the Instructions for Use Lit. 860601. The manufacturer cannot be held liable for instrument cleaning and sterilisation processes carried out at the buyer's facilities.

Important:

For further information about the use of this product please contact your local sales office or visit the homepage www.implan-tec.at. The medical author of the CCG® is available to provide instruction to interested surgeons.

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